

Required Tool: (Not included)

- 4 mm Allen wrench

⚠ WARNING


⚠ CAUTION: To avoid potential shock or other injury, turn power switch **OFF** and disconnect charger before conducting any assembly or maintenance procedures. Failure to follow these steps in the correct order may cause irreparable damage.

Step 1

Using a 4 mm Allen wrench, remove the four (4) bolts from the seat. Remove the seat from the unit.

**Step 2**

Using a 4 mm Allen wrench, remove the four (4) bolts from the seat belt bracket.



Step 3

Remove the seat belt from the unit.



Step 4

Reverse the steps to install the new seat belt.

ATTENTION: Charge battery 12 hours before using.